

Datasheet

HTATIP2 monoclonal antibody (M04), clone 3E11

Catalog Number: H00010553-M04

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a full-length recombinant HTATIP2.

Clone Name: 3E11

Immunogen: HTATIP2 (AAH15358, 1 a.a. ~ 242 a.a.)
full-length recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

MAETEALSKLREDFRMQNKSVFILGASGETGRVLLKEI
LEQGLFSKVTLIGRRKLTDFDEEAYKNVNQEVVDFEKLD
DYASAFQGHVDVGFCCLTTRGKAGAEGFVRVDRDYV
LKSAELAKAGGCKHFNLLSSKGADKSSNFLYLQVKGE
VEAKVEELKFDRYSVFRPGVLLCDRQESRPGEWLVVRK
FFGSLPDSWARGHSVPVTVVRAMLNNVVRPRDKQM
ELLENKAIHDLGKAHGSLKP

Host: Mouse

Reactivity: Human

Applications: ELISA, S-ELISA

(See our web site product page for detailed applications
information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 10553

Gene Symbol: HTATIP2

Gene Alias: CC3, FLJ26963, SDR44U1, TIP30