



Certificate of Conformity

VWR CELL CULTURE FLASK, NON-TREATED

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Cat. No.	Description	Working Volume	Growth Area	Qty/Pk
734-2783	25 ml, Tissue Culture Flask, Non Treated, w/plug seal cap, Sterilized	20 ml	12.5 cm ²	50
734-2784	25 ml, Tissue Culture Flask, Non Treated, w/vent cap, Sterilized	20 ml	12.5 cm ²	50
734-2785	50 ml, Treated for Increased Cell Attachment, w/plug seal cap, Sterilized	40 ml	25.0 cm ²	50
734-2786	50 ml, Treated for Increased Cell Attachment, w/vent cap, Sterilized	40 ml	25.0 cm ²	50
734-2787	250 ml, Treated for Increased Cell Attachment, w/plug seal cap, Sterilized	175 ml	75.0 cm ²	25
734-2788	250 ml, Treated for Increased Cell Attachment, w/vent cap, Sterilized	175 ml	75.0 cm ²	25

Materials:

Flask: GPPS (General polystyrene) Colour: Clear


Cap: HDPE(High Density Polyethylene) Colour: Purple 259

Membrane: 0.22µm Hydrophobic PP (Only use for Vented Cap)








Certificate of Conformity

734-2783 and 734-2784




Dimensions (Unit: mm)






Certificate of Conformity

734-2785 and 734-2786





Dimensions (Unit: mm)





Certificate of Conformity

734-2787 and 734-2788



Other Characteristics:

- Temperature: -20°C~+50°C
- Sterilized by gamma irradiation SAL10⁻⁶(ISO11137)
- Non-pyrogenic
- Bulk packed in zip-Lock bag
- Autoclavability: No
- Every product is printed with a batch number for traceability
- Shelf Life: 3 years after month of production

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **March 2017**