

Certificate of Specification

Material : 84893.290 Kieselgel 60 - 200 µm

CHARACTERISTICS	SPECIFICATIONS
Bulkdichtheid (getast)	460 --> 550 g/l
pH (5 %)	6,5 --> 7,5
Verlies bij drogen (160°C; 15 min)	2,0 --> 4,0 %
Deeltjesgrootte (< 63 µm)	Max. 5,0 %
Deeltjesgrootte (> 200 µm)	Max. 5,0 %

We certify that this batch conforms to the specifications listed above.
BDL : Below detected limit.

Wendy Mathues, Head of laboratory - Haasrode
VWR International
Document printed on 02/2015

This document has been produced electronically and is valid without a signature.